

# INTEGRATING THE HEALTHCARE ENTERPRISE INTEGRATION STATEMENT

## CAPSULE MEDICAL DEVICE INTEGRATION SOLUTION (FORMERLY SMARTLINX)

### Overview

This Integrating the Healthcare Enterprise (IHE) Integration Statement describes the intended conformance of the Capsule product with the IHE Technical Framework. The product implements all transactions required in the IHE Technical Framework to support the IHE Integration and Content Profiles, Actors and Options listed below.

This IHE Integration Statement does not guarantee successful interoperability of the Capsule product with products from other vendors. To ensure safe and effective operation of connected, interoperable products, users of the Capsule product are responsible for testing and validating the full range of product interactions supported by the implemented integration. Testing and validation are strongly recommended by Capsule and other vendors of IHE conforming products, before interoperable products are used clinically.

Capsule accepts that the IHE Technical Framework and the underlying communication standards will evolve and that some advancements will require changes to already installed systems. The company is committed to modifying and adapting its products to maintain conformance to meet market conditions and user requirements.

CAPSULE TECHNOLOGIES IHE INTEGRATION STATEMENT		Date
		2020 - Feb
Vendor	Product Name	Version
Capsule Technologies	Capsule Medical Device Integration Solution (Formerly SmartLinx)	Version v10.4 and v11.0
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:		
Integration Profiles Implemented	Actors Implemented	Options Implemented
Consistent Time	Time Client	None
Device Enterprise Communication	Device Observation Consumer Device Observation Reporter	Waveform Content Module
Alarm Communication Management	Alert Reporter Alert Consumer	Waveform Content Module None
Patient Demographics Query	Patient Data Supplier	None
Links to Standards Conformance Statements for the Implementation		
	<a href="https://capsuletech.com/ihe-statements/ihe-statements-for-product-releases">https://capsuletech.com/ihe-statements/ihe-statements-for-product-releases</a>	
Contact IHE	North America <a href="http://www.ihe.net">www.ihe.net</a> Europe <a href="http://www.ihe-europe.org">www.ihe-europe.org</a> Japan <a href="http://www.ihe-j.org/">http://www.ihe-j.org/</a>	

### Disclaimer

Capsule is not responsible for changes in products from other vendors that affect interoperability. Refer to complete instructions for use documentation for all connected, interoperable products to confirm supported IHE Profiles.

### FOR MORE INFORMATION, CONTACT US

#### NORTH AMERICA

+1 800-260-9537  
support@capsuletech.com

#### INTERNATIONAL OFFICES

+33 1 84 17 12 50  
international@capsuletech.com

Regarding Capsule IHE Conformance, go online to <https://capsuletech.com/ihe-statements/ihe-statements-for-product-releases>