

INTEGRATING THE HEALTHCARE ENTERPRISE CAPSULE MEDICAL DEVICE INTEGRATION SOLUTION INTEGRATION STATEMENT

Overview

This Integrating the Healthcare Enterprise (IHE) Integration Statement describes the intended conformance of the Capsule product with the IHE Technical Framework. The product implements all transactions required in the IHE Technical Framework to support the IHE Integration and Content Profiles, Actors and Options listed below.

This IHE Integration Statement does not guarantee successful interoperability of the Capsule product with products from other vendors. To ensure safe and effective operation of connected, interoperable products, users of the Capsule product are responsible for testing and validating the full range of product interactions supported by the implemented integration. Testing and validation are strongly recommended by Capsule and other vendors of IHE conforming products, before interoperable products are used clinically.

Capsule accepts that the IHE Technical Framework and the underlying communication standards will evolve and that some advancements will require changes to already installed systems. The company is committed to modifying and adapting its products to maintain conformance to meet market conditions and user requirements.

CAPSULE MEDICAL DEVICE INTEGRATION SOLUTION IHE INTEGRATION			
Vendor	Product Name	Version	Date
Capsule Technologies	Capsule Medical Device Integration Solution (formerly SmartLinx)	10.4 and higher	2020 - Feb
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below.			
Profiles Implemented	Actors Implemented	Options Implemented	
Consistent Time	Time Client	None	
Device Enterprise Communication	Device Observation Consumer Device Observation Reporter	Waveform Content Module Waveform Content Module	
Alarm Communication Management	Alert Reporter Alert Consumer	Waveform Content Module None	
Patient Demographics Query	Patient Demographics Supplier	None	
Links to Conformance Statements and Standards			
Capsule	www.capsuletech.com/ihe-integration		
IHE	North America www.ihe.net Europe www.ihe-europe.net Japan www.ihe-j.org		

Disclaimer

Capsule is not responsible for changes in products from other vendors that affect interoperability. Refer to complete instructions for use documentation for all connected, interoperable products to confirm supported IHE Profiles.

FOR MORE INFORMATION, CONTACT US

NORTH AMERICA
+1 800-260-9537
support@capsuletech.com

INTERNATIONAL OFFICES
+33 1 84 17 12 50
international@capsuletech.com

CapsuleTech.com